

Genomic Signature Tags (GST) Protocol

Required reagents and recommended suppliers

T4 DNA ligase 350U/μl from Takara Biotechnology equals ~2.8 Weiss U/μl

T4 DNA ligase High Concentration 5 Weiss U/μl from Invitrogen/Gibco

*Nla*III 10U/μl from New England Bio Labs (NEB) - **Store at -80°C.**

We also obtain *Not*I and *Bam*HI from NEB and
Exonuclease I 20U/μl

*Mme*I 2U/μl is obtained from Centrum Transferu technologii, Ul. Grunwaldzka 529, 80-320
Gdansk, Poland (**sole source**)

Streptavidin coated Dynabeads M-280 from Dynal or equivalent

Dynal MPC magnet or other suitable device.

Oligonucleotides are purchased from Integrated DNA Technologies, Inc. Oligonucleotides are
purified by PAGE or HPLC.

10 x OFA(One Phor All buffer, Amersham, Pharmacia)

GlycoBlue 15 mg/ml (Ambion)

Taq DNA polymerase 10x reaction buffer without MgCl₂-Promega catalogue #M190

Platinum Taq DNA Polymerase High Fidelity (Invitrogen)

Preparation of linker cassettes

Oligonucleotides for the linker cassettes are dissolved in ddH₂O to a concentration of 100 pmoles/μl

Pairs are annealed together in the following standard reaction:

P=phosphorylated ; *= 3' amino modified

36 μl top strand oligonucleotide

36 μl bottom strand oligonucleotide

10 μl 10xOFA

18 μl TEsl (10 mM TrisHCl, pH 8.0; 0.1 m M EDTA-Na₃)

100 μl

95°C 2 min

65°C 10 min

37°C 10 min

RT 20 min--check annealing by electrophoresis on a 10% polyacrylamide gel,
use 20-40 pmol of each separate strand and annealed product,
store at -20°C

Preparation of genomic DNA

DNA is digested with *Bam*HI or other suitable fragmenting enzyme. Typical reaction would be

10 μg DNA in a final vol. of 100 μl NEB#2 buffer plus 1xBSA with 1μl enzyme for 2 hrs at 37°

PC (Phenol/chloroform) extracted and EtOH ppt from 0.3 M NaOAc, pH 6.0, Collect by centrifugation
and resuspend in 34 μl TEsl

Ligation of biotinylated matching enzyme linker cassette (*Bam*HI cassette is shown)

5'-CGAACCCCTTCG

BIOTIN-TGCTTGGGGAAGCCTAGp

To 34 μ l enzyme digested genomic DNA, add

5 μ l 10x T4 DNA ligase buffer (Takara)

8 μ l matching enzyme linker cassette (~50 fold excess, ~288 pmoles)

3 μ l T4 DNA ligase (Takara)

50 μ l incubate at 16°C O/N

PC extract to remove / inactivate ligase, wash PC phase with 50 μ l TEsl, pool aqueous phases then ethanol precipitate. Chill sample at -80°C for 1 hour or O/N at -20°C, then spin 30 min at 4°C , wash with cold 75% EtOH, dry.

1st *Nla*III digestion

Rehydrate pellet in 83 μ l ddH₂O, add:

10 μ l NEB#4

1.0 μ l 100xBSA 10 mg/ml

4 μ l 100mM spermidine(HCl)₃

add 2.0 μ l *Nla*III, *Nla*III is stored at -80°C

100 μ l Final vol. incubate 3 hrs. at 37°C

Bind biotinylated fragments to streptavidin beads.

Remove 100 μ l thoroughly resuspended Dynal M280 streptavidin beads from the stock into a clean 1.5 ml siliconized or low adhesion (Ambion) microcentrifuge tube and place tube on magnet. Remove supernatant and wash beads with 400 μ l 1x B&W buffer (binding and wash buffer 10 mM TrisHCl, pH 8.0, 1 M NaCl, 1 mM EDTA-Na₃)

Resuspend beads in 100 μ l **2x** B&W buffer

Add 100 μ l digest

Incubate at RT for 1 hr. with gentle mixing. **Do NOT vortex**, but make sure beads are fully resuspended.

2nd *Nla*III digestion

Collect beads, carefully remove supernatant and add premixed:

168 μ l ddH₂O

20 μ l 10xNEB #4

2 μ l 100xBSA

8 μ l 100mM spermidine(HCl)₃

2 μ l *Nla*III

PCR amplification of GSTs

primers:

forward is biotinylated and corresponds to a portion of *MmeI* linker's top strand

5'-Biotin-GGATTTGCTGGTCGAGTACA

reverse is biotinylated and corresponds to a portion of degenerate linker's bottom strand

5'-Biotin-TAGTCAGTTGCGACACTAGTGGC

GST PCR cycle

95°C 2 min

95°C 30 sec

58°C 30 sec

72°C 30 sec

72°C 4 min

10°C hold

30 cycles steps 2-4

PCR Reaction: use Promega buffer- Do not use a PCR buffer with ammonium sulfate as it is insoluble in EtOH which then causes problems with later steps.

1 x µl	10 x µl	20 x µl	Stock Conc.	Reagent
18.55	185.5	371		ddH ₂ O
2.5	25	50	10x	Promega Buffer
1.0	10.0	20.0	50 mM	MgSO ₄
0.75	7.5	15.0	10 mM each	dNTPs
1.0	10.0	20.0	10 µM	biotinyl forward primer
1.0	10.0	20.0	10 µM	biotinyl reverse primer
0.1	1.0	2.0		Platinum Taq polymerase mix
0.1	1.0	2.0		cDNA products*
25.0 µl	250 µl	500 µl	(five or ten tubes of 50µl)	

* We initially tried a range of concentrations, from 0.01µl to 1.0µl, but found that this produced the best yield - might need to optimize.

Pool products and then use 200µl in following reaction:

Linear amplification to resolve heteroduplexes (LARHD)

1 x µl	20 x µl	40 x µl	50 x µl	Stk Conc	Reagent
13.65	273	546	682.5		ddH ₂ O
2.5	50	100	125	10x	Promega buffer
1.0	20	40	50	50 mM	MgSO ₄
0.75	15	30	37.5	10 mM each	dNTPs
1.0	20	40	50	10 µM	biotinyl forward primer
1.0	20	40	50	10 µM	biotinyl reverse primer
5.0	100	200	250		1 st round amp'd tags
0.1	2	4.0	5		Platinum Taq polymerase mix
25.0 µl	500µl	1000µl	1250µl	(10, 20 or 25 tubes of 50µl)	

LARHD cycle (only one cycle)

95°C 2.5 min
 58°C 30 sec
 72°C 5 min
 10°C hold

Pool products, saving 10µl for gel analysis if desired. Products should be 94 bp in length

[----- **GST** -----]

5' -GGATTTGCTGGTCGAGTACAACCTAGGCTTAATCCG**CATG** NNNNNNNNNNNNNNNNNNNNTT**CATG**GCGGAGA
 3' -CCTAAACGACCAGCTCATGTTGATCCGAATTAGGCT**GTAC** NNNNNNNNNNNNNNNNNNNA**GTAC**CGCCTCT

CGTCCGCCACTAGTGTGCGCAACTGACTA
 GCAGGCGGTGATCACAGCGTTGACTGAT

Exonuclease I digestion of primers

To the pooled PCR products add:

10µl *E. coli* Exonuclease I, 10 U/µl, for 1000 µl of LARHD products

Incubate at 37°C 60 minutes, remove 5 µl for gel

PC extract to remove enzymes, wash PC with a small amount of TEsl, pool, then precipitate in as many tubes as needed:

1 tube = 270µl sample

30µl 3M NaOAc, pH 6.0

750µl 100% EtOH, Place at -80°C for 30 min.

spin down amplicons in **cold**, redissolve in **cold** 0.3 M NaOAc and then reprecipitate in one tube with 2.5 vol. EtOH, wash with 70% EtOH, dry.

LARHD-2 25 cycles of linear amplification with biotinyl forward primer followed by **one** cycle of

amplification with biotinyl reverse primer.

1 x μ l	40 x μ l	Stk Conc	Reagent
17.65	706		ddH ₂ O
2.5	100	10x	Promega buffer
1.0	40	50 mM	MgSO ₄
0.75	30	10 mM each	dNTPs
1.0	40	10 μ M	biotinyl forward primer
1.0	40		1 st round amp'd tags
0.1	4.0		Platinum Taq polymerase mix
25.0 μ l	1000 μ l	(10 tubes of 100 μ l)	

95°C 2 min
95°C 30 sec
58°C 30 sec
72°C 30 sec
72°C 5 min
10°C hold

25 cycles steps 2-4

add 4 μ l biotinyl B reverse primer and an additional 0.4 μ l of Platinum Taq polymerase mix to each 100 μ l reaction mix, followed by one cycle of denaturation and extension.

95°C 2.5 min
58°C 30 sec
72°C 5 min
10°C hold

Exonuclease I digestion of primers

To the pooled PCR products add:

10 μ l *E. coli* Exonuclease I, for 1000 μ l of LARHD-2 products

Incubate at 37°C 60 minutes, remove 5 μ l for gel

PC extract to remove enzymes, wash PC with a small amount of TEsl, pool, then precipitate in as many tubes as needed:

1 tube = 270 μ l sample
30 μ l 3M NaOAc, pH 6.0
750 μ l 100% EtOH, Place at -80°C for 30 min.

spin down amplicons in **cold**, redissolve in **cold** 0.3 M NaOAc, reprecipitate in one tube with 2.5 vol. EtOH, wash with 70% EtOH, dry.

***Na*III digestion** *Na*III is stored at -80°C

Digestion is performed at 37°C for 4 hrs in 400 µl 1 x NEB #4 buffer plus 1x BSA and 4 mM permidine(HCl)₃

*Nla*III (2 µl) is typically added twice, 2x 2 hr digestions for a total of 4 hrs.

Digestion products should be:

***Mme*I linker arm-40mer**

1 **Biotin**-GGATTTGCTG GTCGAGTACA ACTAGGCTTA ATCCGACATG
CCTAAACGAC CAGCTCATGT TGATCCGAAT TAGGCT

plus GST-23mer

RRRRRR RRRRRRRRRR RTTCATG
GTACYYYYYY YYYYYYYYYY YAA

plus degenerate linker arm -35mer

GCGGAG ACGTCCGCCA CTAGTGTCGC AACTGACTA
GTACCGCCTC TGCAGGCGGT GATCACAGCG TTGACTGAT-**Biotin**

PC extract on **ice**, EtOH ppt from 0.3 M NaOAc plus 2.5 µl GlycoBlue, chill at -80°C,
SPIN IN COLD ROOM (IMPORTANT)

wash pellet in **ice-cold** 70% EtOH, dry

Resuspend sample in 200 µl **cold** TEsl +25 mM NaCl

Bind biotinylated arms to Dynal streptavidin beads:

use 200µl Promega beads, prewashed in 1x B&W Buffer (1 M NaCl- no added BSA or glycogen)

resuspend washed beads in 200µl 2x B&W, add *Nla*III digest, mix at RT for 15-30 min,

Collect beads and save unbound fraction --**THESE ARE THE GSTs.**

Wash with 100µl 1x B&W buffer, pool and ppt with 2.5 vol. EtOH at -80°C, **SPIN IN COLD ROOM (IMPORTANT) wash pellet in ice-cold 70% EtOH, dry**

Self-ligation of cassette tags to form concatemers

Resuspend GST pellet in 12.5µl TEsl on ice, add

1.5µl 10 x T4 DNA ligase buffer (Takara)

1.0µl T4 DNA Ligase (Gibco-High Conc)

Incubate at 16°C 4-6 hours

Add 25 µl TEsl+25 mM NaCl, heat for 2.5 min at 65°C, quench on ice, add 4µl 80% glycerol, mix and apply to a single slot of an 8 slot 0.75% Low Melt agarose minigel bracketed by 100 and 500 bp ladders. Electrophoresis and cut out concatemers.

Note in some cases we cut the gel to remove tags of < 100-250 bp, and reversed the gel's polarity to concentrate the DNA prior to elution.

Tags are purified using GFX Spin columns (Amersham, Pharmacia).

Elute with 180 μ l ddH₂O
add 20 μ l 3M NaOAc, pH 6.0,
2 μ l GlycoBlue
500 μ l EtOH, chill spin, wash with 70% EtOH, dry,

up in 8 μ l TEsl
plus 1 μ l 10x T4 DNA ligase buffer Takara, mix, add
0.5 μ l *Sph*1 cut pZero-BNL, heat at 65°C for 30 sec, quench on ice, add
1 μ l T4 DNA ligase (Takara), mix, incubate at 16°C for several hrs.
dilute to 50 μ l by adding 40 μ l 1X T4 Ligation buffer plus 1 μ l T4 DNA ligase, incubate O/N at 16°C,

PC extract, EtOH ppt from NaOAc + GlycoBlue

Sample up in 10 -15 μ l ddH₂O, Electroporate 5 μ l sample into 50 μ l TOP10 competent cells.

Phenotypically express each sample in 1 ml 2xYT for 1 hr at 37°C with shaking, pool = 3 ml.
Plate 200, 100, 50 and 25 μ l onto prewarmed 2xYT plates + 50 μ g/ml kanamycin.
Add 0.1 vol of 80% glycerol to the remaining cells and store at -80°C.
Incubate plates overnight. A good library should provide 100-200 colonies on the 25 μ l to 50 μ l platings.

GST concatemers in

pZero are sequenced with the m13 FORWARD primer

GSTs should have the following sequence where polarity is indicated by R and Y

5' -RRRRRR RRRRRRRRRR **RTTCATG**-3'
3' -**GTAC**YYYYYY YYYYYYYYYY **YAA**-5'

Concatemers should have the following type GST units which are extracted using software developed at BNL

1 RRRRRR RRRRRRRRRR **RTTCATGAAY** YYYYYYYYYY YYYYYY**CATG** RRRRRRRRRR
GTACYYYYYY YYYYYYYYYY **YAAGTACTTR** RRRRRRRRRR RRRRRR**GTAC** YYYYYYYYYY

57 RRRRRR**RTTC** **ATGAAY**YYYYY YYYYYYYYYY **YCATG**RRRRR RRRRRRRRRR **RRRTTCATGR**
YYYYYY**AAG** **TACTT**RRRRR RRRRRRRRRR **RRGTAC**YYYYY YYYYYYYYYY **YYAAGTACY**

117 RRRRRRRRRR RRRRRR**TTCA** **TGAAY**YYYYY YYYYYYYYYY **YCATG**RRRRR RRRRRRRRRR
YYYYYYYYYY YYYYYY**AAGT** **ACTT**RRRRR RRRRRRRRRR **RGTAC**YYYYY YYYYYYYYYY

177 **RRTTCATGR**R RRRRRRRRRR RRRRR**TTCA**T **G**RRRRRRRRR RRRRRRRR**TT** **CATGAAY**YYY
YYAAGTACYY YYYYYYYYYY YYYYY**AAGTA** **C**YYYYYYYYY YYYYYYY**AA** **GTACTT**RRRR

237 YYYYYYYYYY **YYYCATGAAY** YYYYYYYYYY YYYYYY**CATG** **AA**YYYYYYYYY YYYYYYYYY**C**
RRRRRRRRRR **RRRGTACTTR** RRRRRRRRRR RRRRRR**GTAC** **TT**RRRRRRRRR RRRRRRRR**R**

297 **ATGAAY**YYYYY YYYYYYYYYY **YCATG**
TACTTRRRRR RRRRRRRRRR **RR**

***MmeI* Long SAGE Protocol**

This method uses many of the same steps and reagents as in the GST protocol. The main difference is that it uses cDNA prepared by Reverse Transcription of poly(A)⁺ mRNA as the starting material. The efficiency of this step is increased by capturing the poly(A)⁺ mRNA on oligo (dT)₂₅ magnetic beads directly from a cell lysate and by repeating the first-strand synthesis step several times.

The following new reagents that are needed:

- Reverse Transcriptase RNase H-free SuperScript II
- DEPC treated dH₂O
- E. coli* DNA ligase, DNA polymerase and RNase H
- 1st and 2nd strand buffers.
- Glycogen 20 mg/ml

All of these reagents are obtained from Invitrogen

In addition, a Dynal Dynabeads mRNA Direct kit (catalogue #610.11) is recommended. This kit contains reagents for cell lysis and includes the oligo (dT)₂₅ beads.

SuperRNasin is obtained from Ambion

Thoroughly resuspend and remove 100µl suspended Dynal oligo(dT) beads from the stock into a clean 1.5 ml siliconized or low adhesion Ambion microcentrifuge tube and place tube on magnet. Remove supernatant and wash beads with 400µl lysis /binding buffer from the Dynal Dynabeads mRNA Direct kit.

Mix RNA with 500µl Dynal Lysis/Binding buffer supplemented with 10 µg/ml glycogen.

Collect washed beads, add RNA solution to beads, mix, heat to 60 °C for 5 min, cool to room temp for 10 min on bench with occasional mixing (about 1x per min)

Collect beads, wash 2 x with 400µl Dynal wash buffer A (with LiDS) supplemented with 20 µg/ml glycogen.

Wash 3 x with 400µl Dynal wash buffer B (w/o LiDS) + 20 µg / ml glycogen.

Move beads to new tube after 1st wash.

Wash 2 x with 400µl RT 1st strand buffer + 20 µg/ml glycogen + 2µl SuperRNasin.

First Strand cDNA synthesis

- a) Resuspend beads in 25µl RT 1st strand buffer (+SuperRNasin. + glycogen as above).
- b) Incubate at 42°C 2 minutes
- c) Incubate at 37°C 2 minutes. Scrape down beads off sides if necessary, add a mixed containing:
 - 9.0µl DEPC treated water
 - 1.0µl SuperRNasin
 - 5.0µl 5x 1st strand buffer
 - 2.5µl dNTPs (10 mM each) then add RT b/f adding to beads)
 - 5.0µl DTT
 - 2.5µl Superscript II Reverse Transcriptase

Incubate at 37°C 1 hr with gentle mixing.

Heat to 60°C 3 minutes, incubate at 37° 2 minutes, then add an additional 2 µl RT. Incubate at 37°C an additional 1 hr. Repeat one more time. Total of 6.5 µl RT is used for 3 cycles of cDNA synthesis.

Collect beads, carefully remove 1st stand reagents.

Second Strand cDNA synthesis

Add to beads (premixed)

- 253.5µl dH₂O
- 70µl 5x 2nd strand buffer
- 8µl dNTPs (10 mM each)
- 2.5µl E. coli DNA ligase
- 10µl E. coli DNA polymerase
- 2.5µl E.coli RNase H
- 3.5µl glycogen @ 5 µg /µl
- 350µl

Mix, incubate at 16°C O/N, with occasional mixing for the first 3 hrs.

- a) Wash beads 6 x with 1x B&W buffer plus 1x BSA (0.1 mg/ml) using 500µl/wash
 - 1st time: Resuspend in 1x B&W buffer / BSA and heat to 75°C for 10 minutes, cool
 - 2nd time: Wash without heat
 - 3-6 times: Quick rinse without heat. After 4th wash, transfer to clean tube.
(can store in 500µl 1x B&W buffer+BSA at 4°C)

Digestion with *Nla*III anchoring enzyme

Wash beads 3 x with 200µl 1x NEB #4 buffer

Resuspend in

200µl 1x NEB #4 buffer plus 1x BSA plus 4 mM spermidine(HCl)₃; add

2µl *Nla*III (10 U/µl from NEB) ***Nla*III is stored at -80°C**

Incubate at 37°C for 2 hours.

add an additional 2µl *Nla*III and incubate for an additional 2 hrs with occasional mixing

Capture beads, remove supernatant and EtOH ppt for gel analysis (digested cDNA should produce a visible smear of low MW fragments ~100 -500 bp)

Resuspend beads in 400µl 1x B&W buffer + BSA and heat at 65°C for 20 min to inactivate *Nla*III; cool and wash 4 more times with 400µl of 1x B&W buffer + BSA. EtOH ppt these washes as well for gel analysis

Store beads in 400µl 1x B&W buffer+ BSA at 4°C.

Proceed to

LIGATION of 1ST LINKER CASSETTE (LINKER A) with site for *Mme*I tagging enzyme.

All subsequent steps are the same as in the GST protocol.